Pressure control valve kit Installation instructions

1. Intended use

The pump and the accessory kit are intended exclusively for use as follows:

- For transferring and metering liquids.
- For temporary operation in laboratories for research purposes.
- For operation in accordance with the operating parameters specified in the technical data and other applicable documents.
- For transferring media which meet the requirements in the technical data and other applicable documents.
- For operation in fully assembled condition.
- For upright operation on a firm table top.

2. Reasonably foreseeable misuse

The pump and the accessory kit may not:

- be operated in an explosive atmosphere.
- be used to transfer explosive media.
- be used to transfer media with a flashpoint of less than 40 °C.
- be used to transfer foodstuffs or pharmacological products. Additional certifications are required for this purpose.
- be used to transfer liquid bromine.
- be used to transfer media whose compatibility with the parts coming into contact with liquid has not been verified.
- be used to transfer media containing solids with a particle size exceeding 70 µm.

3. Safety

⚠ DANGER

Danger of injuries and property damage caused by hazardous materials

Poisoning and caustic burns, or undesirable reactions caused by escaping hazardous materials

- Observe the safety data sheets for the media to be transferred
- Clarify the resistance of the parts coming into contact with liquid
- Maintain the tightening torque of connections
- Check tightening torque of head screws regularly
- Check that pump and system are leak-tight
- Do not use the pump and the pressure control valve unless they are in technically perfect condition
- Always operate the pump and the pressure control valve in accordance with the technical data

⚠ WARNING

Danger of undetected diaphragm rupture

The diaphragm may be damaged by aggressive media

- Flush pressure control valve thoroughly after use or if it is not going to be used for a long period
- Pay attention to the service life of the diaphragm

Before operating the pump and the accessories, please read the operating instructions on the web site (www.knf.com/downloads) and pay attention to the safety precautions!

⚠ CAUTION**Danger of cavitation or particles in the fluid circuit**

Deviating metering accuracy


- Use a prefilter, especially for media containing particulate matter or fibres
- Vent the pump and system thoroughly

⚠ CAUTION**Danger of excessive banking-up pressure (resistance) in the filter as a result of high particle concentration in the medium or a mesh size which is too small**


This reduces the flow rate and causes errors in metering accuracy

- Use filters with a sufficiently large mesh size
- Filter fluid beforehand where necessary
- Check filter penetration regularly, cleaning or replacing the filter where necessary (see filter data sheet)


4. Install components


1. Fit adapter plate (1) to underside of pump housing (3) with countersunk screws (2)


2. Remove protective plug from outlet of pump (4) and inlet of pressure control valve (5)

- Screw 70 mm hose section (6) into pump head (7) and tighten finger-tight (tighten a little more if necessary)


- Screw 70 mm hose section (6) into pressure control valve (8) and tighten finger-tight (tighten a little more if necessary)


- Screw pressure control valve (8) to adapter plate (1) with the two small screws (9). Ensure that the flow direction (arrow) of the pressure control valve is correct.

5. Connect UNF1/4"-28

- Push clamp ring (2) onto flanged hose
- Push hose connector (3) onto flanged hose
- Finger-tighten** the hose connector (3).
- Check that the hoses (4) and transition joints are fitted correctly and securely.
- Check system for leaks

6. Screw fitting with cutting ring and union nut

- Slide union nut (4) onto hose
- Slide cutting ring (5) onto hose
- Slide grommet (6) onto hose
- Slide hose (1) onto connecting piece (2) and tighten with union nut (4) - finger-tight
- Check that the hoses (1) and transition joints (hose connector / hose) are connected correctly and securely
- Check that the system is leak-tight

7. Cleaning and storing the system

In order to maintain the service life of the pump and the pressure control valve, after each use and before any long periods of inactivity the pump must be flushed through with neutral media at atmospheric pressure.

IMPORTANT:

To ensure satisfactory start-up when the unit is started up again, it is important to ensure that the pump is free of any crystallising, adhesive or curing media.

- Flush pump thoroughly (see operating instructions)

8. Check system for leaks

⚠ WARNING


Danger of fluid system bursting due to overpressure.

The pump builds up pressure. In a closed system the max. permitted operating pressure may be exceeded, thus causing damage to the pump and/or the system.

- In the case of parts that are in contact with the fluid, only select those that are designed for the pump's operating pressure
- Do not transfer media against shut-off elements or closed systems, otherwise
- install pressure reduction / safety valve and set to 6 bar (available from KNF)

In order to ensure correct assembly and hence reliability, it is essential to perform a leak test.

1. Set up the system as shown above
2. Use distilled water as test medium
3. Check that screw fitting is not leaking
4. Operate pump carefully for 5 minutes
5. Check pump and pressure control valve for leaks

If there are no traces of medium on the pump head, the pump and the pressure control valve can be used in accordance with the specification.

Caution!

If it is found that liquid is escaping, under no circumstances should the pump and the pressure control valve continue to be used.

1. Check that connections are leak-tight.
2. Check tightening torque of head screws.
3. Check that elastomer parts in pump head are clean.
4. Check leak-tightness again.

If these actions have no effect, contact your local KNF dealer and do not continue to use the pump under any circumstances!

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KNF FLODOS AG