



KNF Neuberger Laboratory Pump Kits

KNF Pump/Model No.	Kit No.	Qty/Pump	Kit Description
N920 Vacuum Pumps			
N920 AP.18	057456	1	(3) EPDM diaphragms, 6 O-rings, (6) EPDM valve plates, (2) O-rings
N920 AP.29.18	057456	1	(3) EPDM diaphragms, (6) O-rings, (6) EPDM valve plates, (2) O-rings
N920 KT.29.18	111905	1	(3) FPM diaphragms, (6) O-rings, (6) FFPM valve plates, (2) O-rings
LABOPORT Vacuum Pumps			
N810 FTP	072778	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N820 ANP	073375	1	(1) CR diaphragm, (2) NBR valves plates/seal rings
N820 FTP	072779	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N840 FTP	072780	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N840.1.2 FTP	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N810.3 FTP	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N820.3 ANP	073375	2	(1) CR diaphragm, (2) NBR valves plates/seal rings
N820.3 FTP	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N840.3 FTP	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N842.3 FTP	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N860.1.2 FTP	075097	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N860.3 FTP	075097	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
Mini LABOPORT Vacuum Pumps			
N86 KTP	072839	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N86 KNP	072838	1	(1)- EPDM diaphragm, (2) FPM valves plates/seal rings
N811 KTP	074931	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N811 KNP	074626	1	(1) EPDM diaphragm, (2) FPM valves plates/seal rings
N811 KVP	074178	1	(1) PTFE coated diaphragm, (2) FPM valves plates/seal rings
N816.3 KTP	077769	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N816.1.2 KTP	077769	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N838.3 KTP	202884	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N838.1.2 KTP	202884	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings



KNF Neuberger Laboratory Pump Kits

KNF Pump/Model No.	Kit No.	Qty/Pump	Kit Description
N86 KT.45P	072839	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N86 KT.46P	072839	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N811 KV.45P	074178	1	(1) PTFE coated diaphragm, (2) FPM valves plates/seal rings
N811 KT.46P	074931	1	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N816.3 KT.45P	077769	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N816.1.2 KT.45P	077769	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N838.3 KT.45P	202884	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N838.1.2 KT.45P	202884	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
POWERDRY Self-Cleaning Vacuum Pumps			
N820.3 FT.40P	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N840.3 FT.40P	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N842.3 FT.40P	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
N860.3 FT.40P	075097	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
Programmable Vacuum Systems			
SC920	111905	1	(3) FPM diaphragms, (6) O-rings, (6) FFPM valve plates, (2) O-rings
SC950	126111	1	(4) PTFE coated diaphragms, (14) FFPM valve plates, (2) seal rings, (24) FPM and FFPM O-rings
LABOPORT SC Vacuum Systems			
SC 810	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SC 820	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SC 840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SC 842	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
LABOXACT Manual Vacuum Systems			
SEM810	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SEM820	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SEM840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SEM842	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings



KNF Neuberger Laboratory Pump Kits

KNF Pump/Model No.	Kit No.	Qty/Pump	Kit Description
LABOBASE Vacuum Systems			
SBC840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SBC840.40	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SBC844	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SBC844.40	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SBC860	075097	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SBC860.40	075097	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
LABOPORT SR Vacuum Systems			
SR 810	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SR 820	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SR 840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SR 842	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
LABOPORT SH Vacuum Systems			
SH 810	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SH 820	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SH 840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SH 842	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
LABOPORT SCC Vacuum Systems			
SCC 810	072778	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SCC 820	072779	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SCC 840	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings
SCC 842	072780	2	(1) PTFE coated diaphragm, (2) FFPM valves plates/seal rings



KNF Neuberger Laboratory Pump Kits

KNF Pump/Model No.	Kit No.	Qty/Pump	Kit Description
N035 Vacuum Pumps			
N035 ANP	071738	1	(1) CR diaphragm, (2) SS reed valves, (1) CR head gasket, (2) screws
N035 ATP	071739	1	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035 AVP	071740	1	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035 STP	071734	1	(1) PTFE coated diaphragm, (1) PTFE valve plate, (1) PTFE protection disc
N035 SVP	071735	1	(1) FPM diaphragm, (1) FPM valve plate
N035 TTP	071739	1	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035 TVP	071740	1	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035.1.2 ANP	071738	2	(1) CR diaphragm, (2) SS reed valves, (1) CR head gasket, (2) screws
N035.1.2 ATP	071739	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.1.2 AVP	071740	2	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035.1.2 STP	071734	2	(1) PTFE coated diaphragm, (1) PTFE valve plate, (1) PTFE protection disc
N035.1.2 SVP	071735	2	(1) FPM diaphragm, (1) FPM valve plate
N035.1.2 TTP	071739	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.1.2 TVP	071740	2	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035.3 ANP	071738	2	(1) CR diaphragm, (2) SS reed valves, (1) CR head gasket, (2) screws
N035.3 ATP	071739	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.3 AVP	071740	2	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035.3 STP	071734	2	(1) PTFE coated diaphragm, (1) PTFE valve plate, (1) PTFE protection disc
N035.3 SVP	071735	2	(1) FPM diaphragm, (1) FPM valve plate
N035.3 TTP	071739	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.3 TVP	071740	2	(1) FPM diaphragm, (2) SS reed valves, (1) FPM head gasket, (2) screws
N035 TTP - HI FLOW	071737	1	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.1.2 TTP - HI FLOW	071737	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws
N035.3 TTP - HI FLOW	071737	2	(1) PTFE coated diaphragm, (2) SS reed valves, (1) FPM head gasket, (1) PTFE protection disc, (2) screws



KNF Neuberger Laboratory Pump Kits

KNF Pump/Model No.	Kit No.	Qty/Pump	Kit Description
N022 Vacuum Pumps			
N022 ANP	071723	1	(1) CR diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) CR head gasket, (1) screw
N022 ATP	071724	1	(1) PTFE coated diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) PTFE head gasket, (1) screw
N022 AVP	071725	1	(1) FPM diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) FPM head gasket, (1) screw
N022 STP	071726	1	(1) PTFE coated diaphragm, (2) SS reed valves, (1) PTFE head gasket, (2) screw
N022 SVP	071722	1	(1) FPM coated diaphragm, (1) FPM valve plate
N026 Pump Series			
N026 ANP	071723	1	(1) CR diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) CR head gasket, (1) screw
N026 ATP	071724	1	(1) PTFE coated diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) PTFE head gasket, (1) screw
N026 AVP	071725	1	(1) FPM diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) FPM head gasket, (1) screw
N026 STP	071721	1	(1) PTFE coated diaphragm, (1) PTFE valve plate
N026 SVP	071722	1	(1) FPM coated diaphragm, (1) FPM valve plate
N026.1.2 Pump Series			
N026.1.2 ANP	071723	2	(1) CR diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) CR head gasket, (1) screw
N026.1.2 ATP	071724	2	(1) PTFE coated diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, PTFE head gasket, (1) screw
N026.1.2 AVP	071725	2	(1) FPM diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) FPM head gasket, (1) screw
N026.1.2 STP	071721	2	(1) PTFE coated diaphragm, (1) PTFE valve plate
N026.1.2 SVP	071722	2	(1) FPM coated diaphragm, (1) FPM valve plate
N026.3 Pump Series			
N026.3 ANP	071723	2	(1) CR diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) CR head gasket, (1) screw
N026.3 ATP	071724	2	(1) PTFE coated diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, PTFE head gasket, (1) screw
N026.3 AVP	071725	2	(1) FPM diaphragm, (2) SS reed valves, (2) washers, (1) hex nut, (1) FPM head gasket, (1) screw
N026.3 STP	071721	2	(1) PTFE coated diaphragm, (1) PTFE valve plate
N026.3 SVP	071722	2	(1) FPM coated diaphragm, (1) FPM valve plate